

Product Datasheet

GPR12 Antibody (Extracellular Domain) (orb1533894)

Description	GPR12 Antibody (Extracellular Domain)
Species/Host	Rabbit
Reactivity	Bovine, Equine, Hamster, Human, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 15 amino acid peptide from 2nd extracellular domain of human GPR12 / GPCR12. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Marmoset, Mouse, Rat, Hamster, Elephant, Panda, Bovine, Horse, Rabbit, Pig (100%); Monkey, Dog, Bat, Opossum, Turkey, Chicken, Xenopus (93%); Platypus, Pufferfish (87%); Stickleback (80%).
Target	GPR12
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1.2 µg/ml)
Expiration Date	12 months from date of receipt.