

Product Datasheet

SLC7A5 / CD98 Light Chain Antibody (Internal) (orb1533896)

Description	SLC7A5 / CD98 Light Chain Antibody (Internal)
Species/Host	Rabbit
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 15 amino acid peptide from internal region of human SLC7A5. Percent identity with other species by BLAST analysis: Monkey, Marmoset, Dog (100%); Human, Panda, Bat, Rabbit (93%); Elephant, Bovine, Horse (87%).
Target	SLC7A5 / CD98 Light Chain
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	0.6 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (7 µg/ml)
Expiration Date	12 months from date of receipt.