

Product Datasheet

GRM3 / MGLUR3 Antibody (N-Terminus) (orb1533917)

Description	GRM3 / MGLUR3 Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IHC-P, WB
Immunogen	Synthetic 17 amino acid peptide from N-Terminus extracellular domain of human GRM3 / MGLUR3. Percent identity with other species by BLAST analysis: Human, Gorilla, Orangutan, Gibbon, Monkey, Marmoset, Elephant (100%); Rat, Hamster, Panda (94%); Mouse, Horse, Rabbit, Opossum (88%); Bat (82%).
Target	GRM3 / MGLUR3
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ICC, IHC, IHC-P (9 µg/ml), WB
Expiration Date	12 months from date of receipt.