

Product Datasheet

HTR1F / 5-HT1F Receptor Antibody (N-Terminus) (orb1533937)

Description	HTR1F / 5-HT1F Receptor Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Bovine, Human, Monkey, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 19 amino acid peptide from N-Terminus extracellular domain of human 5HT1F Receptor. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Rat, Bovine, Rabbit (100%); Mouse, Elephant, Panda, Bat, Guinea pig (95%); Platypus (84%).
Target	HTR1F / 5-HT1F Receptor
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 - 10 µg/ml)
Expiration Date	12 months from date of receipt.