

Product Datasheet

HTR1F / 5-HT1F Receptor Antibody (Cytoplasmic Domain) (orb1533939)

Description

HTR1F / 5-HT1F Receptor Antibody

Species/Host	Rabbit
Reactivity	Bat, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 19 amino acid peptide from 3rd cytoplasmic domain of human 5HT1F Receptor. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Monkey, Panda, Bat (100%); Gibbon, Mouse, Rat, Hamster, Elephant, Dog, Bovine, Rabbit (95%); Horse (89%).
Target	HTR1F / 5-HT1F Receptor
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (6 - 9 µg/ml)
Expiration Date	12 months from date of receipt.