

Product Datasheet

GALR1 / Galanin Receptor 1 Antibody (C-Terminus) (orb1534071)

Description	GALR1 / Galanin Receptor 1 Antibody (C-Terminus)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 20 amino acid peptide from C-Terminus of human GALR1. Percent identity with other species by BLAST analysis: Gorilla (100%); Human, Chimpanzee, Orangutan, Gibbon (95%); Monkey (89%).
Target	GALR1 / Galanin Receptor 1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (7 µg/ml)
Expiration Date	12 months from date of receipt.