

Product Datasheet

GALR1 / Galanin Receptor 1 Antibody (N-Terminus) (orb1534096)

Description	GALR1 / Galanin Receptor 1 Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from N-Terminus extracellular domain of human GALR1. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey (100%); Marmoset (94%); Dog (89%); Panda, Horse (83%).
Target	GALR1 / Galanin Receptor 1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ICC, IHC, IHC-P (19 µg/ml)
Expiration Date	12 months from date of receipt.