

Product Datasheet

RORB / ROR Beta Antibody (Hinge Domain) (orb1534278)

Description

RORB / ROR Beta Antibody (Hinge)

Species/Host	Rabbit
Reactivity	Bat, Bovine, Canine, Equine, Gallus, Hamster, Human, Monkey, Mouse, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from hinge domain of human ROR Beta. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Mouse, Rat, Sheep, Hamster, Panda, Bovine, Dog, Bat, Horse, Rabbit, Opossum, Chicken, Platypus, Lizard (100%); Elephant, Xenopus (94%); Zebrafish (83%).
Target	RORB / ROR Beta
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml)
Expiration Date	12 months from date of receipt.