

Product Datasheet

CELSR1 Antibody (Extracellular Domain) (orb1534286)

Description	CELSR1 Antibody (Extracellular Domain)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 15 amino acid peptide from 1st extracellular domain of human CELSR1. Percent identity with other species by BLAST analysis: Human, Gorilla, Elephant, Panda (100%); Gibbon, Mouse, Hamster, Bovine, Dog, Horse (93%); Rabbit (87%); Platypus (80%).
Target	CELSR1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (15 - 20 µg/ml)
Expiration Date	12 months from date of receipt.