

Product Datasheet

HTR1B / 5-HT1B Receptor Antibody (Cytoplasmic Domain) (orb1534312)

Description

HTR1B / 5-HT1B Receptor

Species/Host

Rabbit

Reactivity

Bat, Bovine, Canine, Guinea pig, Hamster, Human, Mouse, Porcine, Rabbit, Rat

Conjugation

Unconjugated

Tested

IHC, IHC-P

Applications
Immunogen

Synthetic 16 amino acid peptide from 3rd cytoplasmic domain of human 5HT1B Receptor. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Marmoset, Mouse, Rat, Hamster, Elephant, Panda, Bat, Bovine, Dog, Cat, Rabbit, Pig, Guinea pig (100%); Monkey, Horse, Opossum (94%); Chicken, Lizard, Xenopus (88%).

Target

HTR1B / 5-HT1B Receptor

Preservatives

PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide

Concentration

1 mg/ml

Storage

Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (7 µg/ml)

Expiration Date

12 months from date of receipt.