


Product Datasheet

HTR7 / 5HT7 Receptor Antibody (N-Terminus) (orb1534345)

Description

HTR7 / 5HT7 Receptor Antibody (N-


Species/Host

Rabbit

Reactivity

Canine, Human, Monkey

Conjugation

Unconjugated

Tested


IHC, IHC-P

Applications

Immunohistochemical
staining of human
br...

Immunogen

Synthetic 20 amino acid peptide from N-Terminus extracellular domain of human 5HT7 Receptor. Percent identity with other species by BLAST analysis: Human, Marmoset, Dog (100%); Mouse, Rat, Hamster, Panda (95%); Gorilla, Pig (90%); Horse (85%).


Target

HTR7 / 5HT7 Receptor

Preservatives

PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide

Immunohistochemical
staining of human
br...

Form/Appearance

PBS, 0.1% Sodium Azide

Concentration

p>0.6 mg/ml

Storage

Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Purity

Immunoaffinity purified

Dilution Range

IHC, IHC-P (10 µg/ml)

Expiration Date

12 months from date of receipt.