

Product Datasheet

ESRRB / ERR Beta Antibody (Internal) (orb1534426)

Description	ESRRB / ERR Beta Antibody (Internal)
Species/Host	Rabbit
Reactivity	Equine, Hamster, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 14 amino acid peptide from internal region of human ESRRB. Percent identity with other species by BLAST analysis: Human, Gorilla, Monkey, Marmoset, Hamster, Panda, Horse (100%); Mouse, Rat, Elephant, Bovine, Bat, Rabbit, Pig (93%); Dog (86%).
Target	ESRRB / ERR Beta
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1 µg/ml)
Expiration Date	12 months from date of receipt.