

Product Datasheet

TRPV4 Antibody (N-Terminus) (orb1534642)

Description	TRPV4 Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from N-Terminus of human TRPV4. Percent identity with other species by BLAST analysis: Human (100%); Marmoset (94%); Panda, Dog (89%); Hamster, Bovine, Pig (83%).
Target	TRPV4
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	0.91 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (2.5 - 5 µg/ml)
Expiration Date	12 months from date of receipt.