

Product Datasheet

SLC12A2 / NKCC1 Antibody (Internal) (orb1534713)

Description	SLC12A2 / NKCC1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from internal region of human SLC12A2 / NKCC1. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset (100%); Bovine, Pig (89%); Ferret, Dog, Bat, Panda (83%).
Target	SLC12A2 / NKCC1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (7.5 - 10 µg/ml)
Expiration Date	12 months from date of receipt.