

Product Datasheet

SLC5A1 / SGLT1 Antibody (N-Terminus) (orb1534736)

Description	SLC5A1 / SGLT1 Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 16 amino acid peptide from N-Terminus of human SLC5A1 / SGLT1. Percent identity with other species by BLAST analysis: Human (100%); Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey (94%); Baboon (88%).
Target	SLC5A1 / SGLT1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 - 15 µg/ml)
Expiration Date	12 months from date of receipt.