

Product Datasheet

DPP4 / CD26 Antibody (Internal) (orb1534740)

Description	DPP4 / CD26 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 20 amino acid peptide from internal region of human CD26. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon (100%); Chimpanzee, Orangutan, Monkey (95%); Galago (90%); Marmoset, Panda, Horse, Rabbit, Pig (85%); Elephant, Cat, Dog (80%).
Target	DPP4 / CD26
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 - 10 µg/ml)
Expiration Date	12 months from date of receipt.