

Product Datasheet

SLC11A2 / DMT1 Antibody (Internal) (orb1534796)

Description	SLC11A2 / DMT1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 16 amino acid peptide from internal region of human SLC11A2. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Marmoset (100%); Galago, Rat, Hamster, Bovine, Horse (94%); Mouse, Elephant (88%); Dog, Bat, Rabbit, Pig, Guinea pig (81%).
Target	SLC11A2 / DMT1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	0.9 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (15 - 20 µg/ml)
Expiration Date	12 months from date of receipt.