

Product Datasheet

ITGA4 / VLA-4 / CD49d Antibody (C-Terminus) (orb1534826)

Description	ITGA4 / VLA-4 / CD49d Antibody (C-Terminus)
Species/Host	Rabbit
Reactivity	Bat, Canine, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from C-Terminus of human ITGA4 / Integrin Alpha 4. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Mouse, Rat, Dog, Bat, Hamster, Panda, Rabbit, Pig (100%); Bovine, Elephant, Horse, Opossum, Platypus (94%).
Target	ITGA4 / VLA-4 / CD49d
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (2.5 - 5 µg/ml)
Expiration Date	12 months from date of receipt.