

Product Datasheet

VR1 / TRPV1 Antibody (Internal) (orb1534925)

Description	VR1 / TRPV1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Guinea pig, Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 17 amino acid peptide from internal region of human TRPV1. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Guinea pig (100%); Monkey (94%); Horse (82%).
Target	VR1 / TRPV1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml)
Expiration Date	12 months from date of receipt.