

Product Datasheet

HLA-A/B/C Antibody (clone W6/32) (orb1535033)

Description

HLA-A/B/C Antibody (clone W6/32)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC, IHC-Fr

Applications
Immunogen

Human MHC Class I HLA-A/HLA-B/HLA-C

Target

HLA-A/B/C

Preservatives

PBS, pH 7.2, 150 mM NaCl, 0.09% Sodium Azide. PBS, pH 7.2, 150 mM NaCl, 0.09% Sodium Azide

Concentration

0.5 mg/ml

Storage

Store at 4°C, avoid freeze-thaw cycles.

Note

For research use only

Application notes

Not recommended for: IHC-P. Further information: The W6/32 antibody has been tested by flow cytometric analysis of human peripheral blood leukocytes. This can be used at less than or equal to 1 ug per test. A test is defined as the amount (ug) of antibody that will stain a cell sample in a final volume of 100 ul. Cell number should be determined empirically but can range from 10⁵ to 10⁸ cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.

Isotype

IgG2a,k

Clonality

Monoclonal

Dilution Range

Flo, IHC, IHC-Fr (10 µg/ml)

Expiration Date

12 months from date of receipt.