

Product Datasheet

NOS2 / iNOS Antibody (clone K13-A) (orb1535074)

Description	NOS2 / iNOS Antibody (clone K13-A)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-Fr, IHC-P
Immunogen	Peptide derived from human iNOS sequence
Target	NOS2 / iNOS
Preservatives	20 mM Tris-HCl, pH 8.0, 20 mg/ml rAlbumin stabilizer, 0.05% NaN ₃ preservative
Storage	Store at 4°C. Avoid freezing.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-Fr, IHC-P (1:200)
Expiration Date	12 months from date of receipt.