

Product Datasheet

KRT20 / CK20 / Cytokeratin 20 Antibody (clone E19-I) (orb1535076)

Description	KRT20 / CK20 / Cytokeratin 20 Antibody (clone E19-I)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Peptide derived from C-Terminus of human cytokeratin 20
Target	KRT20 / CK20 / Cytokeratin 20
Preservatives	20 mM Tris-HCl, pH 8.0, 20 mg/ml rAlbumin stabilizer, 0.05% NaN ₃ preservative
Storage	Store at 4°C. Avoid freezing.
Note	For research use only
Application notes	Further information: Centrifuge the vial before use.
Clonality	Polyclonal
Dilution Range	IHC-Fr, IHC-P (1:100)
Expiration Date	12 months from date of receipt.