

Product Datasheet

GLB1 / Beta-Galactosidase Antibody (clone 5H2) (orb1535080)

Description	GLB1 / Beta-Galactosidase Antibody (clone 5H2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IHC-P, WB
Immunogen	Full length human recombinant protein of human GLB1(NP_001073279) produced in HEK293T cell.
Target	GLB1 / Beta-Galactosidase
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol, 1% rAlbumin
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Flo (1:100), IF (1:100), IHC, IHC-P (10 µg/ml), WB (1:1000)
Expiration Date	12 months from date of receipt.