



Product Datasheet

PAPLN Antibody (orb153510)


Description	Mouse monoclonal antibody to PAPLN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human PAPLN (AA: 766-870) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	137.7kDa
Entrez	89932
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human PAPLN using the antibody.



Western blot analysis of HEK293 cells (Lane 1) using the antibody.