

Product Datasheet

**LAMP1 / CD107a Antibody (aa32-350, clone 5H6)
(orb1535116)**

Description

LAMP1 / CD107a Antibody (aa32-350,

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IHC-P, WB
Immunogen	Recombinat protein encoding aa 32-350 of NP_005552.3 expressed in E. coli.
Target	LAMP1 / CD107a
Preservatives	PBS, 10 mM Sodium Azide. PBS, 10 mM Sodium Azide
Concentration	1 mg/ml
Storage	Store at 4°C or -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The antibody solution can be used at dilutions of at least 1:2000 in immunofluorescence experiments. In western blotting using chemiluminescence it can be used at dilutions of 1:10000 or lower.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	ICC, IF (1:2000), IHC, IHC-P (10 µg/ml), WB (1:10000)
Expiration Date	12 months from date of receipt.