

Product Datasheet

RNF4 Antibody (aa95-123) (orb1535139)

Description	RNF4 Antibody (aa95-123)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	KLH-conjugated synthetic peptide between amino acids 95-123 from C-Terminus of human RNF4.
Target	RNF4
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2°C to 8°C for up to 6 months. For long term storage store at -20°C.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1:100), WB (1:1000)
Expiration Date	12 months from date of receipt.