

Product Datasheet

HOXA9 Antibody (aa245-272) (orb1535202)

Description	HOXA9 Antibody (aa245-272)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P, WB
Target	HOXA9
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Concentration	2 mg/ml
Storage	Maintain refrigerated at 2°C to 8°C for up to 6 months. For long term storage store at -20°C.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Flo (1:10 - 1:50), IHC, IHC-P (1:75), WB (1:1000)
Expiration Date	12 months from date of receipt.