

Product Datasheet

HOMER1 / Homer 1 Antibody (aa72-339) (orb1535207)

Description

HOMER1 / Homer 1 Antibody (aa72-

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IHC-P, WB
Immunogen	Recombinant protein encompassing a sequence within the center region of human Homer1. The exact sequence is proprietary.
Target	HOMER1 / Homer 1
Preservatives	1X PBS, pH 7, 0.025% ProClin 300, 20% Glycerol
Concentration	0.75 mg/ml
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short term storage (1-2 weeks), store at 4°C. For long term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Further information: IHC-paraffin: Suggested antigen retrieval using heat mediated 10 mM Citrate buffer (pH 6.0).
Clonality	Polyclonal
Dilution Range	ICC (1:100 - 1:1000), IHC, IHC-P (7.5 µg/ml), WB (1:500 - 1:3000)
Expiration Date	12 months from date of receipt.