

Product Datasheet

SCG2 / Secretogranin II Antibody (orb1535227)

Description

SCG2 / Secretogranin II Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications
Immunogen

Recombinant protein encompassing a sequence within the N-terminus region of human Chromogranin C. The exact sequence is proprietary.

Target

SCG2 / Secretogranin II

Preservatives

1X PBS, pH 7, 0.01% Thimerosal, 20% Glycerol, 1% rAlbumin

Concentration

1 mg/ml

Storage

Keep as concentrated solution. For short term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Application notes

Further information: Suggested starting dilutions are as follows: Western blot: 1:500-1:3000, IHC-Paraffin: 7.5 µg/ml. Not yet tested in other applications. Optimal working dilutions should be determined experimentally by the end user.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (7.5 µg/ml), WB (1:500 - 1:3000)

Expiration Date

12 months from date of receipt.