

## Product Datasheet

### AIF1 / IBA1 Antibody (orb1535277)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | AIF1 / IBA1 Antibody                                                 |
| <b>Species/Host</b>        | Rabbit                                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                    |
| <b>Conjugation</b>         | Unconjugated                                                         |
| <b>Tested Applications</b> | ELISA, IHC, IHC-P, WB                                                |
| <b>Target</b>              | AIF1 / IBA1                                                          |
| <b>Preservatives</b>       | PBS, pH 7.4, 0.02% Sodium Azide. PBS, pH 7.4, 0.02% Sodium Azide     |
| <b>Concentration</b>       | 1 mg/ml                                                              |
| <b>Storage</b>             | Store at -20°C for up to 1 year.                                     |
| <b>Note</b>                | For research use only                                                |
| <b>Clonality</b>           | Polyclonal                                                           |
| <b>Dilution Range</b>      | ELISA (1:10000 - 1:40000), IHC, IHC-P (5 µg/ml), WB (1:500 - 1:1000) |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                      |