

Product Datasheet

NEFM / NF-M Antibody (aa250-263) (orb1535281)

Full Length Reagents.

NEFM / NF-M Antibody (aa250-263)

Species/Host	Goat
Reactivity	Bovine, Canine, Equine, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Peptide with sequence C-QASHITVERKDYLK, from the internal region of the protein sequence according to NP_005373.2.
Target	NEFM / NF-M
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Peptide ELISA: antibody detection limit dilution 1:128000. Western blot: Approx 170-200kDa band observed in Mouse Brain lysates (calculated MW of 102kDa according to NP_005373.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.03-0.1 ug/ml.
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), Peptide-ELISA (1:128000), WB (0.03 - 0.1 µg/ml)
Expiration Date	12 months from date of receipt.