

Product Datasheet

NAK / TBK1 Antibody (aa514-527) (orb1535309)

Description

NAK / TBK1 Antibody (aa514-527)

Species/Host

Goat

Reactivity

Human, Monkey

Conjugation

Unconjugated

Tested

ELISA, IHC, IHC-P, WB

Applications
Immunogen

Peptide with sequence C-TIETSLQDIDSRLS, from the internal region of the protein sequence according to NP_037386.1.

Target

NAK / TBK1

Preservatives

TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin

Concentration

0.5 mg/ml

Storage

Aliquot and store at -20°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Further information: Peptide ELISA: antibody detection limit dilution 1:128000. Western blot: Approx 85kD band observed in lysates of cell line HeLa (calculated MW of 83.6kD according to NP_037386.1). Recommended concentration: 1-3 ug/ml. An additional band of unknown identity was also consistently observed at 24kD. This band was successfully blocked by incubation with the immunizing peptide.

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (5 µg/ml), Peptide-ELISA (1:128000), WB (1 - 3 µg/ml)

Expiration Date

12 months from date of receipt.