

Product Datasheet

GYPA / CD235a / Glycophorin A Antibody (clone YTH89.1) (orb1535332)

Description	GYPA / CD235a / Glycophorin A Antibody (clone YTH89.1)
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Cytotoxic Assay, FC, IHC, IHC-Fr, IHC-P
Target	GYPA / CD235a / Glycophorin A
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistology: This product does not require antigen retrieval using heat treatment prior to staining of paraffin sections.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	Ctox, Flo (1:10 - 1:25), IHC, IHC-Fr (1:100 - 1:1000), IHC-P (2.5 - 5 µg/ml)
Expiration Date	12 months from date of receipt.