

Product Datasheet

SLC4A3 / AE3 Antibody (aa36-85) (orb1535377)

Description	SLC4A3 / AE3 Antibody (aa36-85)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide located between aa36-85 of human SLC4A3 (F8WD49). Percent identity by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Monkey, Mouse, Rat, Ferret, Dog, Bovine, Bat, Elephant, Horse, Guinea pig (100%); Rabbit (85%).
Target	SLC4A3 / AE3
Preservatives	PBS, 0.09% sodium azide, 2% sucrose. PBS, 0.09% sodium azide, 2% sucrose
Concentration	0.5 mg/ml
Storage	Short term: Store at 2-8°C. Long term: Aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), WB
Expiration Date	12 months from date of receipt.