

Product Datasheet

RAC1 Antibody (C-Terminus) (orb1535383)

Description

RAC1 Antibody (C-Terminus)

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, IHC-P, WB

Applications

Immunogen

Synthetic peptide from C-Terminus of human RAC1 (P60763, NP_008839).
Percent identity by BLAST analysis:
Human, Chimpanzee, Gorilla, Gibbon, Monkey, Galago, Marmoset, Mouse, Rat, Sheep, Elephant, Panda, Dog, Bovine, Bat, Rabbit, Horse, Pig, Opossum, Guinea pig, Turkey, Zebra finch, Chicken, Platypus, Lizard, Xenopus, Trout, Seabass, Catfish, Salmon, Smelt, Stickleback, Pike, Zebrafish (100%); Sea squirt, Drosophila, Slime mold, Nematode, Beetle (92%).

Target

RAC1

Preservatives

PBS, 0.09% sodium azide, 2% sucrose.
PBS, 0.09% sodium azide, 2% sucrose

Concentration

0.5 mg/ml

Storage

Short term: Store at 2-8°C. Long term:
Aliquot and store at -20°C. Avoid
freeze/thaw cycles.

Note

For research use only

Application notes

Further information: Western Blot:
Suggested dilution at 1 ug/ml in 5% skim
milk / PBS buffer, and HRP conjugated
anti-Rabbit IgG should be diluted in 1:
50,000 - 100,000 as secondary antibody.

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (10 µg/ml), WB (0.2 - 1 µg/ml)

Expiration Date

12 months from date of receipt.