

Product Datasheet

CTNND1 / p120 Catenin Antibody (clone EP66) [Out of stock] (orb1535531)

Description	CTNND1 / p120 Catenin Antibody (clone EP66) [Out of stock]
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, IHC-Fr, IHC-P
Target	CTNND1 / p120 Catenin
Preservatives	Cell culture supernatant that is concentrated, dialyzed, filter sterilized and diluted in buffer pH 7.5, containing rAlbumin and 0.02% sodium azide.
Storage	Store at 2°C to 8°C.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Paraffin and frozen. Control tissue: Breast, Breast Lobular Carcinoma. Localization: Membranous, Cytoplasmic.
Isotype	IgG
Clonality	Monoclonal
Dilution Range	IHC, IHC-Fr (1:100 - 1:500), IHC-P (1:100 - 1:500)
Expiration Date	12 months from date of receipt.