

Product Datasheet

AP1S3 Antibody (orb1535536)

Description	AP1S3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Recombinant human AP-1 complex subunit sigma-3 protein.
Target	AP1S3
Preservatives	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.. PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IHC, IHC-P (10 µg/ml), WB
Expiration Date	12 months from date of receipt.