

Product Datasheet

TOMM20 Antibody [Out of stock] (orb1535548)

Description	TOMM20 Antibody [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-145 of human TOMM20 (NP_055580.1). YCIYFDRKRRSDPNFKNRLRERRKKQKLAKERAGLSKLPDLKDAAVQKFFLEEI QLGEELLAQGEYEKGVDHLTNAIAVCGQPQQLLQVLQQTLPVPVFQMLLTCLPTI SQRIVSAQSLAEDDVE
Target	TOMM20
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Concentration	1.74 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 16kDa, while the observed MW by Western blot was 16kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF (1:20 - 1:50), IHC (1:50 - 1:200), IHC-P (1:200), IP (1:20 - 1:50), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.