

Product Datasheet

**IL25 / IL17E Antibody (aa115-132, clone
68C1039.2) (orb1535549)**

Description

IL25 / IL17E Antibody (aa115-132,

Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide from human IL25 / IL-25.
Target	IL25 / IL17E
Preservatives	PBS, 0.05% Sodium Azide. PBS, 0.05% Sodium Azide
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: LS-B123 was validated for use in immunohistochemistry on a panel of 21 formalin-fixed, paraffin-embedded (FFPE) human tissues after heat induced antigen retrieval in pH 6.0 citrate buffer. After incubation with the primary antibody, slides were incubated with biotinylated secondary antibody , followed by alkaline phosphatase-streptavidin and chromogen. The stained slides were evaluated by a pathologist to confirm staining specificity. The optimal working concentration for LS-B123 was determined to be 10 ug/ml.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (2 - 4 µg/ml)
Expiration Date	12 months from date of receipt.