



Product Datasheet

BDNF Antibody (orb153555)


Description	Mouse monoclonal antibody to BDNF
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Purified recombinant fragment of human BDNF (AA: 19-248) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	27.8kDa
Entrez	627
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human BDNF using...



Western blot analysis of HEK293 (Lane 1)...