

Product Datasheet

SLC9A1 / NHE1 Antibody (orb1535600)

Description	SLC9A1 / NHE1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IHC-P, WB
Immunogen	Synthetic peptide derived from human SLC9A1 C-Terminus cytoplasmic region
Target	SLC9A1 / NHE1
Preservatives	PBS, pH 7.2, 50% Glycerol, No preservatives added. PBS, pH 7.2, 50% Glycerol, No preservatives added
Storage	Store at -20°C for up to 1 year.
Note	For research use only
Application notes	Further information: IHC: Heat treatment is necessary for paraffin-embedded sections;
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Flo (1:100), ICC (1:100), IHC, IHC-P (1:100), WB (1:500)
Expiration Date	12 months from date of receipt.