

Product Datasheet

ACE2 / ACE-2 Antibody (aa18-237, clone 1G4) (orb1535606)

Description	ACE2 / ACE-2 Antibody (aa18-237, clone 1G4)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 18-237 of human ACE2 (NP_068576) produced in E. coli.
Target	ACE2 / ACE-2
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol, 1% rAlbumin
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze-thaw cycles. Stable for 12 months from date of receipt if properly stored.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (1:2000)
Expiration Date	12 months from date of receipt.