

Product Datasheet

GC / Vitamin D-Binding Protein Antibody [Out of stock] (orb1535645)

Description	GC / Vitamin D-Binding Protein Antibody [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-493 of human GC (NP_001191236.1). AMDVFCVCTYFMPAAQLPELPDVELPTNKDVCDPGNTKVMDKYTFELSRRTHLPE VFLSKVLEPTLKSGLGECCDVEDSTTCFNAKGPLLKKELSSFIDKQELCADYSENT FTEYKKKLAERLKAKLPDATPTELAKLVNKHSDFASNCCSINSPLYCDSEIDAELK NIL
Target	GC / Vitamin D-Binding Protein
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 39kDa/52kDa/55kDa, while the observed MW by Western blot was 60kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IF (1:10 - 1:100), IHC (1:50 - 1:200), IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.