

Product Datasheet

ANGPTL3 Antibody (aa71-120) (orb1535692)

Description

ANGPTL3 Antibody (aa71-120)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Synthetic peptide located between aa71-120 of human ANGPTL3 (Q9Y5C1, NP_055310). Percent identity by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Galago, Marmoset, Mouse, Panda, Dog, Bovine, Rabbit, Pig (100%); Rat, Elephant, Bat, Horse, Guinea pig (92%); Zebrafish (85%).
Target	ANGPTL3
Preservatives	PBS, 0.09% sodium azide, 2% sucrose. PBS, 0.09% sodium azide, 2% sucrose
Concentration	0.5 mg/ml
Storage	Short term: Store at 2-8°C. Long term: Aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Further information: Western Blot: Suggested dilution at 1µg/ml in 5% skim milk / PBS buffer, and HRP conjugated anti-Rabbit IgG should be diluted in 1:50,000 - 100,000 as second antibody. ELISA titer in peptide based assay: 1:62500.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 µg/ml), WB (0.2 - 1 µg/ml)
Expiration Date	12 months from date of receipt.