

Product Datasheet

OGG1 Antibody (orb1535793)

Description	OGG1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	A synthetic peptide of human OGG1
Target	OGG1
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide. PBS, pH 7.3, 0.02% Sodium Azide
Concentration	1.53 mg/ml
Storage	Store at 4°C, avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 22kDa/36kDa/38kDa/39kDa/40kDa/45kDa/47kDa, while the observed MW by Western blot was 36kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (1:100), WB (1:200 - 1:1000)
Expiration Date	12 months from date of receipt.