

Product Datasheet

GRN / Granulin Antibody [Out of stock] (orb1535808)

Description	GRN / Granulin Antibody [Out of stock]
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-310 of human GRN (NP_002078.1). CPDGQFCPVACCLDPGGASYSCCRPLLDKWPTTLRHLGGPCQVDAHCSAGHS CIFTVSGTSSCCPFPEAVACGDGHHCCPRGFHCSADGRSCFQRSGNNSVGAIQ CPDSQFECPDFSTCCVMVDGSGWCCPMPQASCCEDRVHCCPHGAFCDLVHTR CITPTGTHPLAKKLPAQRTNRAVALSSVMCPDARSRCPDGSTCCELP SGK YGC CPMPNATCCSDHLHCCPQDTVCDLIQSKLSKENATTDLLTKLPAHTVGDVKCD MEVSCPDGYTCCRLQSGAWGCCPFT
Target	GRN / Granulin
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol. PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol
Storage	Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: The predicted MW is 44kDa/46kDa/63kDa, while the observed MW by Western blot was 50kDa.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC (1:50 - 1:200), IHC-P (1:100), WB (1:500 - 1:2000)
Expiration Date	12 months from date of receipt.