

Product Datasheet

**DLK1 / Pref-1 Antibody (Extracellular Domain,
clone PF299-1) (orb1535957)**

Description

DLK1 / Pref-1 Antibody

	1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IHC-P, WB
Immunogen	Recombinant human DLK1 (extracellular domain).
Target	DLK1 / Pref-1
Preservatives	PBS, pH 7.4, no preservatives.. PBS, pH 7.4, no preservatives.
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: ELISA - Direct or indirect: 1:2000-1:10000. Flow Cytometry - 1:100-1:500. Immunohistochemistry - Paraffin sections 10ug/mL. Western blot - 1:2000-1:10000 using ECL. Suggested blocking and dilution buffer is TBST with 0.05% Tween 20 and 5% skim milk. Suggested incubation time is 1 hour at room temperature.
Isotype	IgG1,k
Clonality	Monoclonal
Dilution Range	ELISA (1:2000 - 1:10000), Flo (1:100 - 1:500), IHC, IHC-P (10 µg/ml), WB (1:2000 - 1:10000)
Expiration Date	12 months from date of receipt.