

Product Datasheet

ANGPTL3 Antibody (Coiled-coil domain) (orb1535960)

Description	ANGPTL3 Antibody (Coiled-coil domain)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Recombinant human ANGPTL3.
Target	ANGPTL3
Preservatives	PBS, pH 7.4, no preservatives.. PBS, pH 7.4, no preservatives.
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: ELISA - Direct and indirect: 1:2000-1:5000. Western Blot - 1:2000-1:5000 using ECL. suggested blocking and dilution buffer is PBST containing 0.05% Tween 20 and 5% skim milk. suggested incubation time is 1 hour at room temperature.
Clonality	Polyclonal
Dilution Range	ELISA (1:2000 - 1:5000), IHC, IHC-P (5 µg/ml), WB (1:2000 - 1:5000)
Expiration Date	12 months from date of receipt.