

Product Datasheet

ANGPTL3 Antibody (Fibrogen-like domain) (orb1535961)

Description

ANGPTL3 Antibody (Fibrogen-like)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Fibrinogen-like domain of recombinant human ANGPTL3.
Target	ANGPTL3
Preservatives	PBS, pH 7.4, no preservatives.. PBS, pH 7.4, no preservatives.
Concentration	1 mg/ml
Storage	Short term: store at 4°C. Long term: store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: ELISA - Direct and indirect: 1:2000-1:5000. Western Blot - 1:2000-1:5000 using ECL. suggested blocking and dilution buffer is PBST containing 0.05% Tween 20 and 5% skim milk. suggested incubation time is 1 hour at room temperature.
Clonality	Polyclonal
Dilution Range	ELISA (1:2000 - 1:5000), IHC, IHC-P (5 µg/ml), WB (1:2000 - 1:5000)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558