

Product Datasheet

LYPLA2 Antibody (orb1535993)

Description	LYPLA2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IHC-P
Immunogen	Recombinant human Acyl-protein thioesterase 2 protein.
Target	LYPLA2
Preservatives	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.. PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Concentration	1 mg/ml
Storage	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IF, IHC, IHC-P (10 µg/ml)
Expiration Date	12 months from date of receipt.